

Product Datasheet

c-Kit Antibody / CD117 / SCFR (orb1823455)

Description This mAb recognizes a protein of 145kDa, identified as CD117/p145kit. It is found on a wide variety ...

Species/Host Mouse

Reactivity Human

Conjugation Unconjugated

Tested Applications IHC-P

Immunogen Recombinant full-length human protein was used as the immunogen for the CD117 antibody.

Storage Aliquot the CD117 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.

Note For research use only

Formula 0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide

Isotype Mouse IgG1, kappa


Clonality Monoclonal

Uniprot ID **P10721**


Hazard Information This CD117 antibody is available for research use only.

Dilution Range Immunohistochemistry (FFPE): 2-4ug/ml for 30 min at RT

Expiration Date 12 months from date of receipt.



Analysis of a HuProt (TM) microarray con...



IHC staining of FFPE human GIST with CD1...