



Product Datasheet

S100P Antibody / MIG9 (orb1823466)


Description	S100P is a 95-amino-acid protein and a member of the S100 family. S100P has been shown to mediate tu...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	A recombinant partial protein sequence (within amino acids 1-95) from the human protein was used as the immunogen for the S100P antibody.
Storage	Aliquot the S100P antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG2b, kappa
Clonality	Monoclonal
Uniprot ID	P25815
Hazard Information	This S100P antibody is available for research use only.
Dilution Range	Western blot: 1-2ug/ml, Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.



Western blot testing of human SK-BR-3 ce...



IHC staining of FFPE human lung cancer w...



IHC staining of FFPE human placental tis...