


Product Datasheet

Human CD93/C1qR ELISA Kit (orb182347)

Description	Human CD93/C1qR ELISA Kit
Reactivity	Human
Range	78 pg/ml - 5,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: A24-K580
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:100, concentration in serum and plasma is 72-267 ng/ml
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA), saliva and urine.
MW	68560 MW
Uniprot ID	Q9NPY3
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human CD93/C1qR PicoKine ELISA Kit stand...