



Product Datasheet

PGP9.5 / UchL1 Antibody (orb1823505)

Description	This mAb reacts with a protein of 20-30kDa, identified as PGP9.5, also known as ubiquitin carboxyl-t...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant full-length human protein was used as the immunogen for the PGP9.5 antibody.
Storage	Aliquot the PGP9.5 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Rabbit IgG, kappa
Clonality	Recombinant
Uniprot ID	P09936
Hazard Information	This PGP9.5 antibody is available for research use only.
Dilution Range	Western blot: 2-4ug/ml, Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human cerebellum ti...



SDS-PAGE analysis of purified, BSA-free ...