



Product Datasheet

p63 Antibody / Tumor protein 63 (orb1823523)


Description	p63 is a homolog of the tumor suppressor p53. It is identified in basal cells in the epithelial laye...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A recombinant fragment corresponding to the C-terminus of human TP63 protein was used as the immunogen for the TP63 antibody.
Storage	Aliquot the TP63 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Rabbit IgG, kappa
Clonality	Recombinant
Uniprot ID	Q9H3D4
Hazard Information	This TP63 antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human skin tissue w...



IHC staining of FFPE human prostate canc...



IHC staining of FFPE human prostate canc...