


Product Datasheet

Mouse GPNMB/Osteoactivin ELISA Kit (orb182354)

Description	Mouse GPNMB/Osteoactivin ELISA Kit
Reactivity	Mouse
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: K23-N502
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	63676 MW
Uniprot ID	Q99P91
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Mouse
Osteoactivin/GPNMB
PicoKine ELISA ...