



Product Datasheet

FOXL2 Antibody (orb1823588)


Description	FOXL2 is a critical factor essential for ovary differentiation and maintenance, and repression of th...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Recombinant full-length human FOXL2 protein was used as the immunogen for the FOXL2 antibody.
Storage	Aliquot the FOXL2 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	1 mg/ml in 1X PBS; BSA free, sodium azide free
Isotype	Rabbit IgG, kappa
Clonality	Recombinant
Uniprot ID	P58012
Hazard Information	This FOXL2 antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human parathyroid t...



IHC staining of FFPE human ovarian carci...



SDS-PAGE analysis of purified, BSA-free ...