



Product Datasheet

ZXDC Antibody / ZXD-like zinc finger protein (orb1823650)


Description	ZXDC cooperates with CIITA to promote transcription of MHC class I and MHC class II genes [UniProt].
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FACS, IF, IHC-P
Immunogen	Recombinant full-length human ZXDC protein was used as the immunogen for the ZXDC antibody.
Storage	Aliquot the ZXDC antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG1
Clonality	Monoclonal
Uniprot ID	Q2QGD7
Hazard Information	This ZXDC antibody is available for research use only.
Dilution Range	Flow cytometry: 1-2ug/million cells, Immunofluorescence: 1-2ug/ml, Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.



Immunofluorescent staining of PFA-fixed ...



IHC staining of FFPE human prostate tiss...



Flow cytometry testing of PFA-fixed huma...