

Product Datasheet

PDCD1 Antibody / PD-1 / PD1 (orb1823702)

Description PDCD-1 (programmed cell death-1 protein), also designated CD279, is a type I transmembrane receptor ...

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications ELISA, IHC-P

Immunogen A recombinant fragment corresponding to the extracellular domain of human PDCD1 protein was used as the immunogen for the PD1 antibody.

Storage Aliquot the PD1 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.

Note For research use only

Formula 1 mg/ml in 1X PBS; BSA free, sodium azide free

Isotype Rabbit IgG, kappa


Clonality Recombinant

Uniprot ID [Q15116](#)


Hazard Information This PD1 antibody is available for research use only.

Dilution Range ELISA (Order BSA-free format for coating), Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT

Expiration Date 12 months from date of receipt.



IHC staining of FFPE human tonsil tissue...



IHC staining of FFPE human tonsil tissue...