


## Product Datasheet

### INHA Antibody / Inhibin alpha (orb1823741)

<b>Description</b>	It recognizes a 47kDa protein, which is identified as alpha sub-unit of Inhibin. It is a gonadal pro...
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IHC-P
<b>Immunogen</b>	A recombinant fragment of human Inhibin alpha protein (within amino acids 233-362) was used as the immunogen for the Inhibin alpha antibody.
<b>Storage</b>	Aliquot the Inhibin alpha antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Formula</b>	1 mg/ml in 1X PBS; BSA free, sodium azide free
<b>Isotype</b>	Rabbit IgG, kappa
<b>Clonality</b>	Recombinant
<b>Uniprot ID</b>	<b>P05111</b>
<b>Hazard Information</b>	This Inhibin alpha antibody is available for research use only.
<b>Dilution Range</b>	Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
<b>Expiration Date</b>	12 months from date of receipt.



IHC staining of FFPE human testicular carcinoma tissue.



Negative control: IHC testing of FFPE human testicular carcinoma tissue.