
Product Datasheet

Mucin Gastric MUC5AC Antibody (Monoclonal, 45M1) (orb18238)

Description

Mucin Gastric MUC5AC Antibody (Monoclonal,


Species/Host

Mouse

Reactivity

Human, Rabbit, Rat

Conjugation

Unconjugated

Tested Applications

ICC, IHC, WB

Immunogen

Mucin from human ovarian cyst fluid.

Form/Appearance

Lyophilized

Concentration

Adding 1 ml of PBS buffer will yield a concentration of 100 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 1ml of PBS buffer will yield a concentration of 100ug/ml.

Isotype

Mouse IgG1

Clonality


Monoclonal

MW


585570 MW

Uniprot ID
P98088
Expiration Date


12 months from date of receipt.



IHC analysis of Mucin Gastric using anti...



IHC analysis of Mucin Gastric using anti...



IHC analysis of Mucin Gastric using anti...