
Product Datasheet

**IL1RA Antibody / Interleukin-1 Receptor
Antagonist / IL1RN (orb1823877)**

Description

Two forms of interleukin-1, designated IL-1a and

by ...

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested Applications

IHC-P

Immunogen

A recombinant partial protein sequence (within amino acids 1-200) from the human protein was used as the immunogen for the IL1RN antibody.

Storage

Aliquot the IL1RN antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.

Note

For research use only

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide

Isotype

Mouse IgG2, kappa

Clonality

Monoclonal

Uniprot ID

P18510

Hazard Information


This IL1RN antibody is available for research use only.

Dilution Range


Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT

Expiration Date


12 months from date of receipt.



IHC staining of FFPE human tonsil tissue...



IHC staining of FFPE human adrenal gland...



Analysis of a HuProt (TM) microarray con...