

Product Datasheet


**IL1RA Antibody / Interleukin-1 Receptor
Antagonist / IL1RN (orb1823879)**

Description

Two forms of interleukin-1, designated IL-1a and



by ...


Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested Applications

IHC-P


IHC staining of FFPE human tonsil tissue...

Immunogen

A recombinant partial protein sequence (within amino acids 1-200) from the human protein was used as the immunogen for the Interleukin-1 Receptor Antagonist antibody.

Storage

Aliquot the Interleukin-1 Receptor Antagonist antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.


Note

For research use only

Formula

1 mg/ml in 1X PBS; BSA free, sodium azide free

Isotype

Mouse IgG2, kappa

Clonality

Monoclonal

Uniprot ID

P18510


IHC staining of FFPE human adrenal gland...

Hazard Information

This Interleukin-1 Receptor Antagonist antibody is available for research use only.

Dilution Range

Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT


Expiration Date

12 months from date of receipt.

Analysis of a HuProt (TM) microarray con...