



Product Datasheet

TPO Antibody / Thyroid Peroxidase (orb1823899)


Description	Thyropoxidase (TPO) is a 933 amino acid, type I transmembrane glycoprotein that plays a key role i...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A recombinant fragment of human TPO (within amino acids 600-900) was used as the immunogen for the Thyroid Peroxidase antibody.
Storage	Aliquot the Thyroid Peroxidase antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Rabbit IgG, kappa
Clonality	Recombinant
Uniprot ID	P07202
Hazard Information	This Thyroid Peroxidase antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human thyroid carci...



IHC staining of FFPE human thyroid carci...



SDS-PAGE analysis of purified, BSA-free ...