



Product Datasheet

TPO Antibody / Thyroid Peroxidase (orb1823900)

Description	Thyropoxidase (TPO) is a 933-amino acid, type I transmembrane glycoprotein that plays a key role i...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A recombinant fragment of human protein was used as the immunogen for the recombinant TPO antibody.
Storage	Aliquot the recombinant TPO antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	1 mg/ml in 1X PBS; BSA free, sodium azide free
Isotype	Mouse IgG1, kappa
Clonality	Recombinant
Uniprot ID	P07202
Hazard Information	This recombinant TPO antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human thyroid carci...



SDS-PAGE analysis of purified, BSA-free ...