



Product Datasheet

MYCN Antibody (orb1823904)


Description	The v-Myc oncogene, initially identified in the MC29 avian retrovirus, causes myelocytomas, carcinom...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FACS
Immunogen	A recombinant human MYCN protein (within amino acids 335-463) was used as the immunogen for the MYCN antibody.
Storage	Aliquot the MYCN antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG1
Clonality	Monoclonal
Uniprot ID	P04198
Hazard Information	This MYCN antibody is available for research use only.
Dilution Range	ELISA (Order BSA-free format for coating), Flow cytometry: 1-2ug/million cells
Expiration Date	12 months from date of receipt.



Flow cytometry testing of PFA-fixed huma...



Analysis of a HuProt (TM) microarray con...



SDS-PAGE analysis of purified, BSA-free ...