



Product Datasheet

FABP3 Antibody (orb1823979)

Description	The intracellular fatty acid-binding proteins (FABPs) belongs to a multigene family. FABPs are divid...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A recombinant partial protein sequence (within amino acids 1-127) from the human protein was used as the immunogen for the Cardiac Fatty Acid Binding Protein antibody.
Storage	Aliquot the Cardiac Fatty Acid Binding Protein antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Rabbit IgG, kappa
Clonality	Recombinant
Uniprot ID	P05413
Hazard Information	This Cardiac Fatty Acid Binding Protein antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human heart tissue ...



SDS-PAGE analysis of purified, BSA-free ...