



Product Datasheet

PRAME Antibody / OIP-4 (orb1824044)

Description	Functions as a transcriptional repressor, inhibiting the signaling of retinoic acid through the reti...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A recombinant partial protein sequence (within amino acids 100-509) from the human protein was used as the immunogen for the OIP-4 antibody.
Storage	Aliquot the OIP-4 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	1 mg/ml in 1X PBS; BSA free, sodium azide free
Isotype	Rabbit IgG, kappa
Clonality	Recombinant
Uniprot ID	P78395
Hazard Information	This OIP-4 antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human melanoma tiss...



IHC staining of FFPE human testis tissue...