
Product Datasheet

**MX1 Antibody / MxA / Myxovirus resistance protein
1 (orb1824066)**

Description


The Dynamin family of microtubule-associated



micro...

Species/Host

Mouse


Reactivity

Human

Conjugation

Unconjugated

Tested Applications

IHC-P


IHC staining of FFPE human kidney tissue...

Immunogen

A recombinant partial protein sequence (within amino acids 400-592) from the human protein was used as the immunogen for the Myxovirus resistance protein 1 antibody.

Storage

Aliquot the Myxovirus resistance protein 1 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.


Note

For research use only

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide

IHC staining of FFPE human ovarian cancer...

Isotype

Mouse IgG2, kappa

Clonality

Monoclonal

Uniprot ID


P20591

Hazard Information

This Myxovirus resistance protein 1 antibody is available for research use only.

Dilution Range

Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT



SDS-PAGE analysis of purified, BSA-free ...

Expiration Date

12 months from date of receipt.