


## Product Datasheet

# Chromogranin B Antibody / Secretogranin I (orb1824133)

<b>Description</b>	Chromogranins (secretogranins) are acidic glycoproteins that localize within secretory granules of e...
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	A recombinant partial protein sequence (within amino acids 250-450) from the human protein was used as the immunogen for the Secretogranin I antibody.
<b>Storage</b>	Aliquot the Secretogranin I antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Formula</b>	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
<b>Isotype</b>	Mouse IgG, kappa
<b>Clonality</b>	Monoclonal
<b>Uniprot ID</b>	<b>P05060</b>
<b>Hazard Information</b>	This CHGB antibody is available for research use only.
<b>Dilution Range</b>	ELISA (Order BSA-free format for coating)
<b>Expiration Date</b>	12 months from date of receipt.



Analysis of a HuProt (TM) microarray con...