



Product Datasheet

BATF3 Antibody / SNFT (orb1824174)


Description	SNFT, also known as BATF3 (basic leucine zipper transcription factor, ATF-like 3), JUNDM1 or JDP1, i...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FACS, IHC-P
Immunogen	A recombinant partial protein sequence (within amino acids 1-123) from the human protein was used as the immunogen for the SNFT antibody.
Storage	Aliquot the SNFT antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG1, kappa
Clonality	Monoclonal
Uniprot ID	Q9NR55
Hazard Information	This SNFT antibody is available for research use only.
Dilution Range	Flow cytometry: 1-2ug/million cells, Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.



Analysis of a HuProt (TM) microarray con...



Flow cytometry testing of PFA-fixed huma...



SDS-PAGE analysis of purified, BSA-free ...