



Product Datasheet

CD35 Antibody / CR1 (orb1824210)

Description	CD35, also called CR1 and complement receptor 1, is a cell membrane-bound, monomeric glycoprotein on...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A recombinant partial protein sequence (within amino acids 600-900) from the human protein was used as the immunogen for the recombinant CR1 antibody.
Storage	Aliquot the recombinant CR1 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG1, kappa
Clonality	Recombinant
Uniprot ID	P17927
Hazard Information	This recombinant CR1 antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.



IHC staining
of FFPE
human tonsil
tissue...



SDS-PAGE
analysis of
purified,
BSA-free ...