



Product Datasheet

CD34 Antibody (orb1824219)


Description	This mAb recognizes a single chain, transmembrane, heavily glycosylated protein of 90-120kDa, which ...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A recombinant partial protein sequence (within amino acids 100-300) from the human protein was used as the immunogen for the CD34 antibody.
Storage	Aliquot the CD34 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG1, lambda
Clonality	Recombinant
Uniprot ID	P28906
Hazard Information	This CD34 antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human colon cancer ...



IHC staining of FFPE human stomach tissue...



IHC staining of FFPE human uterus tissue...