



Product Datasheet

LRH1 Antibody / NR5A2 (orb1824239)


Description	NR5A2 (nuclear receptor subfamily 5, group A, member 2), also designated B1-binding factor (B1F or B...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FACS, IHC-P
Immunogen	Recombinant human NR5A2 protein was used as the immunogen for the NR5A2 antibody.
Storage	Aliquot the NR5A2 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	1 mg/ml in 1X PBS; BSA free, sodium azide free
Isotype	Mouse IgG2b
Clonality	Monoclonal
Uniprot ID	000482
Hazard Information	This NR5A2 antibody is available for research use only.
Dilution Range	Flow cytometry: 1-2ug/million cells, Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human prostate cancer...



Flow cytometry testing of PFA-fixed huma...



Analysis of a HuProt (TM) microarray con...