



## Product Datasheet

### CD21 Antibody (orb1824263)

|                            |                                                                                                                       |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | CD21 also known as complement C3d receptor, Epstein-Barr virus receptor, and complement receptor typ...               |
| <b>Species/Host</b>        | Rabbit                                                                                                                |
| <b>Reactivity</b>          | Human                                                                                                                 |
| <b>Conjugation</b>         | Unconjugated                                                                                                          |
| <b>Tested Applications</b> | IHC-P                                                                                                                 |
| <b>Immunogen</b>           | A recombinant fragment corresponding to the C-terminal of human CD21 was used as the immunogen for the CD21 antibody. |
| <b>Storage</b>             | Aliquot the CD21 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.                     |
| <b>Note</b>                | For research use only                                                                                                 |
| <b>Formula</b>             | 1 mg/ml in 1X PBS; BSA free, sodium azide free                                                                        |
| <b>Isotype</b>             | Rabbit IgG, kappa                                                                                                     |
| <b>Clonality</b>           | Recombinant                                                                                                           |
| <b>Uniprot ID</b>          | <b>P20023</b>                                                                                                         |
| <b>Hazard Information</b>  | This CD21 antibody is available for research use only.                                                                |
| <b>Dilution Range</b>      | Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT                                                                |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                       |



IHC staining of FFPE human tonsil tissue...



IHC staining of FFPE human tonsil tissue...