
Product Datasheet

CD45RA Antibody (Leukocyte marker) (orb1824275)

Description

CD45R, also designated CD45 and PTPRC, has been


Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested Applications

IHC-P

IHC staining of FFPE human lymph node ti...

Immunogen

Recombinant full-length human PTPRC protein was used as the immunogen for the CD45RA antibody.

Storage

Aliquot the CD45RA antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.

Note

For research use only

Formula

1 mg/ml in 1X PBS; BSA free, sodium azide free

Isotype

Rabbit IgG, kappa

Clonality

Recombinant

Uniprot ID


P08575

Hazard Information

This CD45RA antibody is available for research use only.

Dilution Range


Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT



IHC staining of FFPE human tonsil tissue...

Expiration Date

12 months from date of receipt.



Negative control: IHC testing of FFPE hu...