



Product Datasheet

CD45 Antibody (Leukocyte marker) (orb1824288)

Description	CD45R, also designated CD45 and PTPRC, has been identified as a transmembrane glycoprotein, broadly ...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant full-length human PTPRC protein was used as the immunogen for the CD45 antibody.
Storage	Aliquot the CD45 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG1, kappa
Clonality	Recombinant
Uniprot ID	P08575
Hazard Information	This CD45 antibody is available for research use only.
Dilution Range	Western blot: 1-2ug/ml, Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.



Western blot testing of human spleen tis...



IHC staining of FFPE human spleen tissue...