



Product Datasheet

COX2 Antibody / PTGS2 (orb1824299)

Description	Prostaglandins are a diverse group of autocrine and paracrine hormones that mediate many cellular an...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A recombinant partial protein sequence (within amino acids 400-604) from the human protein was used as the immunogen for the Cyclooxygenase 2 antibody.
Storage	Aliquot the Cyclooxygenase 2 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	1 mg/ml in 1X PBS; BSA free, sodium azide free
Isotype	Rabbit IgG, kappa
Clonality	Recombinant
Uniprot ID	P35354
Hazard Information	This Cyclooxygenase 2 antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human prostate tiss...



IHC staining of FFPE human prostate tiss...