



Product Datasheet

COX2 Antibody / PTGS2 (orb1824304)

Description	Prostaglandins are a diverse group of autocrine and paracrine hormones that mediate many cellular an...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A recombinant human COX2 protein fragment (within amino acids 504-604) was used as the immunogen for the COX2 antibody.
Storage	Aliquot the COX2 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG1, kappa
Clonality	Monoclonal
Uniprot ID	P35354
Hazard Information	This COX2 antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human renal cell ca...



SDS-PAGE analysis of purified, BSA-free ...