



Product Datasheet

CD45RA Antibody (Leukocyte marker) (orb1824305)

Description	CD45R, also designated CD45 and PTPRC, has been identified as a transmembrane glycoprotein, broadly ...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Recombinant full-length human PTPRC protein was used as the immunogen for the CD45RA antibody.
Storage	Aliquot the CD45RA antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	1 mg/ml in 1X PBS; BSA free, sodium azide free
Isotype	Rabbit IgG, kappa
Clonality	Recombinant
Uniprot ID	P08575
Hazard Information	This CD45RA antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human tonsil tissue...



IHC staining of FFPE human tonsil tissue...