

---

## Product Datasheet

### CD45RA Antibody (Leukocyte marker) (orb1824308)

**Description**

CD45R, also designated CD45 and PTPRC, has been


**Species/Host**

Rabbit

**Reactivity**

Human

**Conjugation**

Unconjugated

**Tested Applications**

IHC-P

IHC staining of FFPE human tonsil tissue...

**Immunogen**

Recombinant full-length human PTPRC protein was used as the immunogen for the CD45RA antibody.

**Storage**

Aliquot the CD45RA antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.

**Note**

For research use only

**Formula**

1 mg/ml in 1X PBS; BSA free, sodium azide free

**Isotype**

Rabbit IgG, kappa

**Clonality**

Recombinant

**Uniprot ID**


**P08575**

**Hazard Information**

This CD45RA antibody is available for research use only.

**Dilution Range**


Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT



IHC staining of FFPE human tonsil tissue...

**Expiration Date**

12 months from date of receipt.



Negative control: IHC testing of FFPE hu...