
Product Datasheet

CD45RA Antibody (Leukocyte marker) (orb1824309)

Description

CD45R, also designated CD45 and PTPRC, has been


Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested Applications

IHC-P

IHC staining of FFPE human tonsil tissue...

Immunogen

Recombinant full-length human PTPRC protein was used as the immunogen for the CD45RA antibody.

Storage


Aliquot the CD45RA antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.

Note

For research use only

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide


Isotype

Rabbit IgG, kappa

IHC staining of FFPE human tonsil tissue...

Clonality

Recombinant

Uniprot ID


P08575

Hazard Information

This CD45RA antibody is available for research use only.

Dilution Range

Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT



Negative control: IHC testing of FFPE hu...

Expiration Date

12 months from date of receipt.