



Product Datasheet

MPZ Antibody / Myelin Protein Zero (orb1824338)

Description	Zero, also known as myelin protein zero (MPZ) is a Type 1 integral membrane glycoprotein that mediat...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Recombinant full-length human MPZ protein was used as the immunogen for the MPZ antibody.
Storage	Aliquot the MPZ antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	1 mg/ml in 1X PBS; BSA free, sodium azide free
Isotype	Mouse IgG2a, kappa
Clonality	Monoclonal
Uniprot ID	P25189
Hazard Information	This MPZ antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human colon tissue



SDS-PAGE analysis of purified, BSA-free ...