



Product Datasheet

RXRG Antibody (orb1824343)


Description	Two families of retinoid receptors, RARs and RXRs, have been identified. Retinoic acid receptors (RA...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FACS
Immunogen	Recombinant human RXR gamma protein was used as the immunogen for the RXRG antibody.
Storage	Aliquot the RXR gamma antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG2b
Clonality	Monoclonal
Uniprot ID	P78347
Hazard Information	This RXR gamma antibody is available for research use only.
Dilution Range	Flow cytometry: 1-2ug/million cells
Expiration Date	12 months from date of receipt.



Analysis of a HuProt (TM) microarray con...



Flow cytometry testing of PFA-fixed huma...



SDS-PAGE analysis of purified, BSA-free ...