



Product Datasheet

MPZ Antibody / Myelin Protein Zero (orb1824344)

Description	Zero, also known as myelin protein zero (MPZ) is a Type 1 integral membrane glycoprotein that mediat...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Recombinant full-length human MPZ protein was used as the immunogen for the Myelin Protein Zero antibody.
Storage	Aliquot the Myelin Protein Zero antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	1 mg/ml in 1X PBS; BSA free, sodium azide free
Isotype	Mouse IgG2a, kappa
Clonality	Monoclonal
Uniprot ID	P25189
Hazard Information	This Myelin Protein Zero antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human colon carcino...



SDS-PAGE analysis of purified, BSA-free ...