

---

## Product Datasheet

### HDGF Antibody / Hepatoma Derived Growth Factor (orb1824353)

**Description**

This gene encodes a nuclear protein containing a

oxi...

**Species/Host**

Mouse

**Reactivity**

Human

**Conjugation**

Unconjugated

**Tested Applications**

FACS, IHC-P

**Immunogen**

Recombinant full-length human Hepatoma Derived Growth Factor protein was used as the immunogen for the HDGF antibody.

**Storage**

Aliquot the Hepatoma Derived Growth Factor antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.

**Note**

For research use only

**Formula**

1 mg/ml in 1X PBS; BSA free, sodium azide free

**Isotype**

Mouse IgG2a, kappa

**Clonality**

Monoclonal

**Uniprot ID**

**P51858**

**Hazard Information**


This Hepatoma Derived Growth Factor antibody is available for research use only.

**Dilution Range**


Flow cytometry: 1-2ug/million cells, Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT

**Expiration Date**


12 months from date of receipt.



IHC staining of FFPE human colon carcinoma...



Flow cytometry testing of PFA-fixed huma...



Analysis of a HuProt (TM) microarray con...