

Product Datasheet

S100A2 Antibody (orb1824425)

Description	S100A2 belongs to the S100 family of proteins containing 2 EF-hand calcium-binding motifs. S100A2 ma...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant full-length human protein was used as the immunogen for the S100 Calcium Binding Protein A2 antibody.
Storage	Aliquot the S100°Calcium Binding Protein A2 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	1 mg/ml in 1X PBS; BSA free, sodium azide free
Isotype	Mouse IgG2b
Clonality	Monoclonal
Uniprot ID	P29034
Hazard Information	This S100 Calcium Binding Protein A2 antibody is available for research use only.
Dilution Range	ELISA (Order BSA-free format for coating)
Expiration Date	12 months from date of receipt.

