



Product Datasheet

Leptin Receptor Antibody / LEPR (orb1824599)


Description	Although there is substantial evidence that body weight is physiologically regulated, the molecular ...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A recombinant partial protein sequence (within amino acids 335-425) from the human protein was used as the immunogen for the Ob-R antibody.
Storage	Aliquot the Ob-R antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	1 mg/ml in 1X PBS; BSA free, sodium azide free
Isotype	Mouse IgG2b, kappa
Clonality	Monoclonal
Uniprot ID	P48357
Hazard Information	This Ob-R antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.



IHC staining of
FFPE human
kidney tissue...



IHC staining of
FFPE human
hepatocellula...



SDS-PAGE
analysis of
purified, BSA-
free ...