



Product Datasheet

GSTM3 Antibody / Glutathione S-Transferase mu 3 (orb1824615)

Description	Members of the glutathione S-transferase (GST) family of proteins function in the detoxification of ...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Recombinant full-length human GSTM3 protein was used as the immunogen for the recombinant GSTM3 antibody.
Storage	Aliquot the recombinant GSTM3 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG2b, kappa
Clonality	Recombinant
Uniprot ID	P21266
Hazard Information	This recombinant GSTM3 antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human kidney tissue...



SDS-PAGE analysis of purified, BSA-free ...