


Product Datasheet

Uroplakin 3B Antibody / UPK3B (orb1824705)

Description	The protein encoded by this gene is a member of the transmembrane 4 superfamily, also known as the t...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P
Immunogen	A recombinant partial protein sequence (within amino acids 100-300) from the human protein was used as the immunogen for the Uroplakin 3B antibody.
Storage	Aliquot the Uroplakin 3B antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Rabbit IgG, kappa
Clonality	Recombinant
Uniprot ID	Q9BT76
Hazard Information	This Uroplakin 3B antibody is available for research use only.
Dilution Range	ELISA (Order BSA-free format for coating),Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human bladder tissue...



IHC staining of FFPE human bladder tissue...