



Product Datasheet

Uroplakin 3B Antibody / UPK3B (orb1824707)

Description	The protein encoded by this gene is a member of the transmembrane 4 superfamily, also known as the t...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P
Immunogen	A recombinant partial protein sequence (within amino acids 100-300) from the human protein was used as the immunogen for the Uroplakin 3B antibody.
Storage	Aliquot the Uroplakin 3B antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	1 mg/ml in 1X PBS; BSA free, sodium azide free
Isotype	Rabbit IgG, kappa
Clonality	Recombinant
Uniprot ID	Q9BT76
Hazard Information	This Uroplakin 3B antibody is available for research use only.
Dilution Range	ELISA (Order BSA-free format for coating), Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human bladder tissue...



SDS-PAGE analysis of purified, BSA-free ...