



Product Datasheet

SBSN Antibody / Suprabasin (orb1824722)


Description	Suprabasin (SBSN) is a 247 amino acid secreted protein that is upregulated in differentiating kerati...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A recombinant fragment of the human protein was used as the immunogen for the Suprabasin antibody.
Storage	Aliquot the Suprabasin antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	1 mg/ml in 1X PBS; BSA free, sodium azide free
Isotype	Mouse IgG
Clonality	Monoclonal
Uniprot ID	Q6UWP8
Hazard Information	This Suprabasin antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human skin tissue w...



IHC staining of FFPE human esophagus tis...



IHC staining of FFPE human skin tissue w...