



Product Datasheet

Uroplakin 3B Antibody / UPK3B (orb1824726)


Description	The protein encoded by this gene is a member of the transmembrane 4 superfamily, also known as the t...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P
Immunogen	A recombinant partial protein sequence (within amino acids 100-300) from the human protein was used as the immunogen for the UPK3B antibody.
Storage	Aliquot the UPK3B antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG
Clonality	Monoclonal
Uniprot ID	Q9BT76
Hazard Information	This UPK3B antibody is available for research use only.
Dilution Range	ELISA (Order BSA-free format for coating), Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human bladder tissue...



IHC staining of FFPE human bladder tissue...



Analysis of a HuProt (TM) microarray con...