



Product Datasheet

SBSN Antibody / Suprabasin (orb1824738)


Description	Suprabasin (SBSN) is a 247 amino acid secreted protein that is upregulated in differentiating kerati...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A recombinant fragment of the human protein was used as the immunogen for the Suprabasin antibody.
Storage	Aliquot the Suprabasin antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG1, kappa
Clonality	Monoclonal
Uniprot ID	Q6UWP8
Hazard Information	This Suprabasin antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human skin tissue w...



IHC staining of FFPE human esophagus tis...



SDS-PAGE analysis of purified, BSA-free ...