



Product Datasheet

ICAM1 Antibody / CD54 (orb1824747)

Description	ICAM-1/CD54 is a single chain glycoprotein of Ig supergene family, present on unstimulated endotheli...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A recombinant fragment of the human protein was used as the immunogen for the ICAM-1 antibody.
Storage	Aliquot the ICAM-1 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG2a, kappa
Clonality	Monoclonal
Uniprot ID	P05362
Hazard Information	This ICAM-1 antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human tonsil tissue...



IHC staining of FFPE human tonsil tissue...