



Product Datasheet

ICAM1 Antibody / CD54 (orb1824763)

Description	ICAM-1 is a single chain glycoprotein of Ig supergene family, present on unstimulated endothelial ce...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A recombinant fragment corresponding to the intracellular portion of the human protein was used as the immunogen for the recombinant CD54 antibody.
Storage	Aliquot the recombinant CD54 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Rabbit IgG, kappa
Clonality	Recombinant
Uniprot ID	P05362
Hazard Information	This recombinant CD54 antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human tonsil tissue...



SDS-PAGE analysis of purified, BSA-free ...