


Product Datasheet

CD147 / Emmprin / Basigin Antibody (orb1824801)


Description	Basigin is expressed more intensely on thymocytes than on mature peripheral blood T cells. It is imp...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant full-length human protein was used as the immunogen for the Basigin antibody.
Storage	Aliquot the Basigin antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG
Clonality	Monoclonal
Uniprot ID	P35613
Hazard Information	This Basigin antibody is available for research use only.
Dilution Range	Western blot: 1-2ug/ml, Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.



Analysis of a HuProt (TM) microarray con...



Western blot testing of human THP-1 cell...



IHC staining of FFPE human placental tis...