


## Product Datasheet

# CD147 / Emmprin / Basigin Antibody (orb1824802)


<b>Description</b>	Basigin is expressed more intensely on thymocytes than on mature peripheral blood T cells. It is imp...
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IHC-P, WB
<b>Immunogen</b>	Recombinant full-length human protein was used as the immunogen for the Basigin antibody.
<b>Storage</b>	Aliquot the Basigin antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Formula</b>	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
<b>Isotype</b>	Mouse IgG
<b>Clonality</b>	Monoclonal
<b>Uniprot ID</b>	<b>P35613</b>
<b>Hazard Information</b>	This Basigin antibody is available for research use only.
<b>Dilution Range</b>	Western blot: 1-2ug/ml, Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
<b>Expiration Date</b>	12 months from date of receipt.



Analysis of a HuProt (TM) microarray con...



Western blot testing of human THP-1 cell...



IHC staining of FFPE human placental tis...