



Product Datasheet

SMAD4 Antibody (orb1824854)

| | |
|----------------------------|--|
| Description | Signaling from the ligand-activated membrane receptor serine/threonine kinases to nuclear targets is... |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | IHC-P |
| Immunogen | A recombinant partial protein sequence (within amino acids 400-552) from the human protein was used as the immunogen for the recombinant SMAD4 antibody. |
| Storage | Aliquot the recombinant SMAD4 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles. |
| Note | For research use only |
| Formula | 1 mg/ml in 1X PBS; BSA free, sodium azide free |
| Isotype | Rabbit IgG, kappa |
| Clonality | Recombinant |
| Uniprot ID | Q13485 |
| Hazard Information | This recombinant SMAD4 antibody is available for research use only. |
| Dilution Range | Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT |
| Expiration Date | 12 months from date of receipt. |



IHC staining of FFPE human pancreatic carcinoma...



SDS-PAGE analysis of purified, BSA-free ...