



## Product Datasheet

### CD147 / Emmprin / Basigin Antibody (orb1824857)


<b>Description</b>	BSG is expressed more intensely on thymocytes than on mature peripheral blood T cells. It is importa...
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, IHC-P
<b>Immunogen</b>	Recombinant full-length human protein was used as the immunogen for the CD147 antibody.
<b>Storage</b>	Aliquot the CD147 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Formula</b>	1 mg/ml in 1X PBS; BSA free, sodium azide free
<b>Isotype</b>	Mouse IgG
<b>Clonality</b>	Monoclonal
<b>Uniprot ID</b>	<b>P35613</b>
<b>Hazard Information</b>	This CD147 antibody is available for research use only.
<b>Dilution Range</b>	ELISA (Order BSA-free format for coating), Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
<b>Expiration Date</b>	12 months from date of receipt.



IHC staining of FFPE human colon carcino...



IHC staining of FFPE human placental tis...



IHC staining of FFPE human renal cell ca...